


Laryngoscopes



1/3


LS-03-414
MILLER
 Laryngoscope Set, Consists Of Battery Handle (Without Battery)
 With Five Miller Laryngoscope Blades, Fig 1 to 5



LS-03-414 A Fig. 0 Working Length, 55 mm (Baby's)
LS-03-414 B Fig. 1 Working Length, 80 mm (Children's)
LS-03-414 C Fig. 2 Working Length, 130 mm (Adolescent's)
LS-03-414 D Fig. 3 Working Length, 170 mm (Women's)
LS-03-414 E Fig. 4 Working Length, 180 mm (Men's)



LS-03-415
GUEDEL-NEGUS
 Laryngoscope Set, Consists Of Battery Handle (Without Battery)
 With Four Miller Laryngoscope Blades



LS-03-415 A Fig. 1 Working Length, 70 mm (Children's)
LS-03-415 B Fig. 2 Working Length, 90 mm (Adolescent's)
LS-03-415 C Fig. 3 Working Length, 110 mm (Women's)
LS-03-415 D Fig. 4 Working Length, 135 mm (Men's)



Laryngoscopes



LS-03-416
 McINTOSH
 Laryngoscope Set, Consists Of Battery Handle (Without Battery)
 With Three McIntosh Blades, Fig. 1, Fig. 2, Fig. 3

1/3




Fig. 0



Fig. 1




Fig. 2




Fig. 3





Fig. 4

LS-03-417
 McINTOSH
 Laryngoscope Set, Consists Of Battery Handle (Without Battery)
 With Five McIntosh Blades, Fig. 0, Fig. 1, Fig. 2, Fig. 3, Fig. 4

LS-03-417 A Fig. 0 Working Length, 55 mm (Baby's)
LS-03-417 B Fig. 1 Working Length, 75 mm (Children's)
LS-03-417 C Fig. 2 Working Length, 90 mm (Adolescent's)
LS-03-417 D Fig. 3 Working Length, 110 mm (Women's)
LS-03-417 E Fig. 4 Working Length, 135 mm (Men's)

LS-03-418
 FOREGGER
 Laryngoscope Set, Consists Of
 Battery Handle (Without Battery)
 With Three Foregger Blades,
 Fig. 1, Fig. 2, Fig. 3,



1/3



Fig. 1



Fig. 2



Fig. 3



Fig. 4

LS-03-419
 FOREGGER
 Laryngoscope Set, Consists Of
 Battery Handle (Without Battery)
 With Four Foregger Blades,
 Fig. 1, Fig. 2, Fig. 3, Fig. 4,

LS-03-419 A Fig. 1 Working Length, 70 mm (Children's)
LS-03-419 B Fig. 2 Working Length, 90 mm (Adolescent's)
LS-03-419 C Fig. 3 Working Length, 110 mm (Women's)
LS-03-419 D Fig. 4 Working Length, 135 mm (Men's)



LS-03-420
 Fiber Optic
 Battery Handle
 30mm diam.



LS-03-421
 Fiber Optic
 Battery Handle
 18mm diam.



LS-03-422
 Replacement Lamp
 (Krypton) 2.5 V



LS-03-423
 Replacement Lamp
 (Halogen) 2.5 V

